



Tongue Scraper

Fresh Breath, Simple Routine. · SKU KD-OC-TS

Tongue scraper series for removing bacteria and debris from the tongue surface to reduce bad breath — available in several handle designs with custom colours and packaging.

Tongue scrapers complete our oral-care range. A smooth scraping edge removes the coating of bacteria and food debris from the tongue that contributes to bad breath, and the comfortable handle makes daily use easy.

SUPER DEER



Specifications

Series	Tongue scraper series
Customisation	Handle colour, packaging and brand printing to order

Features

Smooth scraping edge · Comfortable handle · Custom colours & packs · Clean-room production · Bulk or retail packing

How to use

- Place the scraper at the back of the tongue and draw it forward with light pressure.
- Rinse and repeat two or three times.
- Wash after use.

Applications

Dental clinics, pharmacies, retail distributors and private-label brands.

Packaging & shipping

Series	Packing
Tongue scraper	Individual blister or bulk, to order

Compliance

ISO 13485, CE. Manufactured under an ISO 13485 quality management system. Certificates and test reports available on request.

Request a quotation: <https://kangdemed.com/request-quote?product=tongue-scraper> · Jason.wangxiaokang@gmail.com · WhatsApp +86 199 5663 5178