



Interdental Brushes

Precision Cleaning Where Floss Can't Reach. · SKU KD-OC-IB

I-shaped interdental brushes (from 0.7 mm), L-shaped brushes for back teeth, soft silicone interdental brushes and dual-use toothpick brushes — for cleaning around braces, implants and bridges.

Interdental brushes clean the spaces between teeth and around orthodontic appliances and implants more effectively than floss. The I-shaped brush is the standard straight design, the L-shaped brush angles the head to reach molars, and the silicone brush offers a wire-free, extra-gentle option. The dual-use toothpick brush combines a pick and a brush in one handle.

SUPER DEER



Specifications

Types	I-shaped interdental brush (0.7 mm and up) · L-shaped interdental brush · Silicone interdental brush · Dual-use toothpick brush
Customisation	Wire size, handle colour, packaging and brand printing to order

Features

Sizes from 0.7 mm · I, L and silicone types · Soft, gum-friendly · Private label available · Clean-room production

How to use

- Choose a brush size that fits the interdental space without forcing.
- Insert gently at the gum line and move in and out a few times.
- Rinse after use; replace when the bristles wear.

Applications

Dental clinics, orthodontic practices, pharmacies and private-label brands.

Packaging & shipping

Type	Description
I-shaped interdental brush	Straight head, sizes from 0.7 mm
L-shaped interdental brush	Angled head for back teeth
Silicone interdental brush	Wire-free soft silicone bristles
Dual-use toothpick brush	Pick on one end, brush on the other

Compliance

ISO 13485, CE. Manufactured under an ISO 13485 quality management system. Certificates and test reports available on request.

Request a quotation: <https://kangdemed.com/request-quote?product=interdental-brushes> · Jason.wangxiaokang@gmail.com · WhatsApp +86 199 5663 5178